USER MANUAL

Hereby, MOB, declares that item MO6827 is in Compliance with the essential requirements and other relevant conditions of Directive 2014/30/EU. The full text of the EU declaration of conformity is available at the following internet address: www.momanual.com
MOB, PO BOX 644, 6710 BP (NL).

PO: 41-XXXXXX
Made in China

MO6827

EN

Muscle massage device

Product Specifications

Input Voltage: 5V

Rated Power: 30W

Interface Type: Type C

Main Material: Acrylonitrile butadiene styrene (ABS) Plastic

Use Environment: -5°C-40°C

Battery Capacity: 1200mAh

Power Supply Mode: Rechargeable Lithium Battery

Dimensions: 132mm*134mm*50mm

Product Weight: 550g

Operating Instructions:

A: Switch Indicator

B: Switch Button

C: Charging Port

D: Massager Gun Head

E: Massager Gun Body

1. Before using for the first time ensure the product is fully charged, with a minimum of the three hours charging.
2. Plug in the power cord to the charging port C, make sure to connect the type -C end of the cable to its respective adaptor
3. The switch indicator A will display flashing LED light to indicate that the battery is charging
4. Once the battery is fully charged , the led lights in the switch indicator A will stop flashing
5. Once you are done charging the device , press the switch button B to turn the device on
6. Simply replace the default massage gun head D with your chosen massager head. There are four types of massager heads:

Tips head for all body parts

Small Flat Heads for joints

U-Head for large muscle groups

Ball Head for large muscle groups

Tip Head
for all
body parts

Small Flat
Head
for joints

U-Head
for vertebrae

Ball Head
for large muscle
groups

Warning

1. Do not drop, break, hit or use the product improperly.
2. Do not submerge the product in water
3. Keep devices away from excessive heat, direct sunlight and fire.
4. Avoid excessive force when using the device.
5. The product is for adult use only.

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